

CHALLENGES OF USING STANDARD GASTROSTOMY TUBES FOR LONG-TERM FEEDING



Since their development about 25 years ago, percutaneous gastrostomy tubes have become well-accepted means of providing enteral access for patients with inadequate oral intake.¹

- Widely accepted method of enteral nutrition because of its low morbidity and high success rate.¹

Frequent complications associated with traditional gastrostomy tubes



Peritubal leakage¹



Catheter occlusion¹



Dislodgement¹



Bulky design, posing an adverse cosmetic effect, particularly in the mobile patient²

These complications have paved the way for the development of new approaches to optimal gastrostomy tube type.¹

- Developments such as use of gastrostomy buttons (also called low-profile or skin-level G tubes) were introduced to achieve problem-free long-term enteral feeding.^{1,3}

A 90 to 100% success rate has been shown for placement of gastrostomy buttons.¹

Structure of a typical low-profile gastrostomy tube

- Visible portion of the tube is rod shaped with a cap closure.³
- The tube is premeasured to be flush to the abdomen.³
- Identified by both a French size (outer diameter) and shaft length (measured in cm).³
- Internal retention device may be a water balloon or a non-balloon mushroom-shaped dome.³
- Anti-reflux valve to prevent leakage of gastric contents when uncapping the device.³



Low-profile balloon gastrostomy tube³



Low-profile non-balloon gastrostomy tube³

Advantages of low-profile G tubes over the standard length G tubes

Standard feeding tube vs. low-profile feeding tube³



Standard feeding tube



Low-Profile Feeding Tube




Patients favor its low profile design and more aesthetic nature (suitable for both pediatric and adult patients)^{1,3}



Poses fewer restrictions on mobility³



Does not require tape to secure excess tubing to the skin³



Reduced risk of gastric leakage through the tube (if the cap should open) due to presence of an anti-reflux valve³



Less likely to dislodge by getting caught on something (in case of dislodgement, new buttons can be inserted at the patient's bedside without the need for a return visit to the hospital)^{1,3}



Serves as a useful alternative for adult patients with neurologic disease¹



Buttons are less prone to clogging due to short tube length¹

DID YOU KNOW ?

British Association for Parenteral and Enteral Nutrition (BAPEN) guidelines 2003 state that in patients where cosmetic considerations are important, low-profile 'button' PEGs can be used which contain a built-in anti-reflux valve to prevent leaks when feeding extension tubes are disconnected.⁴

British Society of Gastroenterology (BSG) guidelines 2010 state that in mobile patients, PEG tube may be replaced by a low profile ('button') replacement when the tract is fully developed. Such low-profile devices have also been employed in patients who are at high risk of inadvertent tube displacement.²



OUR SOLUTION

Enteral feeding tubes can make patients feel limited in their lifestyle choices. AVANOS* MIC-KEY* feeding tubes can provide your patients the "cap it and go" option. It offers patients a truly low-profile alternative, a huge "no-show" difference compared to traditional feeding tubes that can show through the clothing, snag, or be pulled out accidentally.⁵



MIC-KEY* G-tube

- Inflatable silicone internal retention balloon⁶
- High performing anti-reflux valve⁶
- SECUR-LOK* extension set connector mechanism with ENFit® Feeding port⁶

References:

1. O'Dowd M, Given MF, Lee MJ. New approaches to percutaneous gastrostomy. Semin Intervent Radiol. 2004;21(3):191-197. 2. Westaby D, Young A, O'Toole P, Smith G, Sanders DS. The provision of a percutaneously placed enteral tube feeding service. Gut. 2010; 59(12):1592-605. 3. Lord LM. Enteral Access Devices: Types, Function, Care, and Challenges. Nutr Clin Pract. 2018; 33(1):16-38. 4. Stroud M, Duncan H, Nightingale J. Guidelines for enteral feeding in adult hospital patients. Gut. 2003; 52(suppl 7):viii-2. 5. Avanos website. Why MIC-KEY* low profile feeding tubes [Internet]. [cited 2021 May 28]. Available from: <https://avanosmedicaldevices.com/digestive-health/patient-experience-digestive-health/why-mic-key-feeding-tubes/>. 6. Product data sheet, MIC* and MIC-KEY* enteral feeding product catalogue 2020