


# ENDOTRACHEAL TUBE POSITIONING: A CRITICAL ISSUE FOR PEDIATRIC AIRWAY MANAGEMENT



Neonates, infants and small children have short tracheas, making endotracheal tube placement difficult.<sup>1</sup>

Correct endotracheal tube positioning in the pediatric patient involves:<sup>2</sup>

- Placement of endotracheal tube (ETT) into the trachea (avoiding endobronchial intubation)<sup>2</sup>
- Achieving proper depth of ETT insertion<sup>2</sup>

## COMPLICATIONS OF MALPOSITIONING OF ETT:

Malpositioning of ETTs within the trachea in pediatric patients is common, reaching a peak of over **35%**.<sup>2</sup>

### Too deep placement<sup>3</sup>

- Injury to carina
- Endobronchial intubation

### Not enough placement depth<sup>3</sup>

- Accidental extubation
- Laryngeal injury, that could cause inadequate ventilation, hypoxia, brain damage and death


## WAYS USED BY CLINICIANS TO CONFIRM PROPER ETT POSITIONING:

- Presence of auscultated breath sounds in all lung fields<sup>2</sup>
- Continuous presence of a normal end-tidal CO<sub>2</sub> (ETCO<sub>2</sub>) curve<sup>2</sup>
- Bilateral chest rise<sup>2</sup>

But, the above ways do not ensure ideal ETT positioning. These indications may be present even when the ETT is too shallow or too deep<sup>2</sup>

## OUR SOLUTION:

The Avanos Pediatric MICROCUFF\* is designed for pediatric anatomy.




Short, cylindrical cuff placed near the tracheal tube tip, secures cuff placement in the trachea and not the pressure-sensitive larynx



Anatomically-based intubation depth mark results in correct placement and a cuff-free subglottic zone

Tube Size I.D	Age/Weight Years/kg
3.0 mm	term/ ≥ 3kg – 8 months
3.5 mm	8 months- 2 years
4.0 mm	2-4 years
4.5 mm	4-6 years
5.0 mm	6-8 years
5.5 mm	8-10 years
6.0 mm	10-12 years
6.5 mm	12-14 years
7.0 mm	14-16 years

Size selection chart allows easy, accurate tube selection, proven 98.4% accurate as cited in a 500 patient study<sup>4</sup>



Four precision bands to facilitate and confirm optimal tube placement

## DID YOU KNOW

- Many of the currently available cuffed pediatric tracheal tubes have intubation depth marks that were absent or misleading with up to three marks (in contrast to Avanos Pediatric MICROCUFF\* that has four precision bands).<sup>5</sup>
- Intubation depth markings allow instant appropriate placement of the tube in children and are superior to the conventional age based formula for oral tube insertion depth.<sup>1</sup>

**References:** 1. Weiss M, Balmer C, Dullenkopf A, Knirsch W, Gerber AC, Bauersfeld U, Berger F. Intubation depth markings allow an improved positioning of endotracheal tubes in children. *Can J Anaesth.* 2005; 52(7):721. 2. Harris EA, Arheart KL, Penning DH. Endotracheal tube malposition within the pediatric population: a common event despite clinical evidence of correct placement. *Can J Anaesth.* 2008; 55(10):685-90. 3. Neunhoeffer F, Wahl T, Hofbeck M, et al. A new method for determining the insertion depth of tracheal tubes in children: a pilot study. *Br J Anaesth.* 2016; 116(3):393-7. 4. Dullenkopf A, Gerber AC, Weiss M. Fit and seal characteristics of a new paediatric tracheal tube with high volume-low pressure polyurethane cuff. *Acta Anaesthesiol Scand.* 2005; 49(2):232-7. 5. Weiss M, Dullenkopf A. Cuffed tracheal tubes in children: past, present and future. *Expert review of medical devices.* 2007; 4(1):73-82.